

Preliminary Studies of Using Matrix-Assisted Laser Desorption/Ionization Time-of-Flight for Detection of *CFTR* Variants

Victoria M. Pratt Ph.D.¹, Andrew Bradford Ph.D.¹, Katarzyna Gawlik Ph.D.¹, Laura Mendez B.S.¹, Priscilla Hunt B.S.¹, Parastou Sadatmousavi M.S.¹, Vineet K. Dhiman Ph.D.², Dylan Hendricks M.S.², Gajendra Katara Ph.D.³, Christine Banaszak M.S.³, Matthias Jost Ph.D.¹, Glenn Sawyer B.S., M.B.A.¹

¹Agena Bioscience, San Diego CA | ²Illinois Department of Public Health, Chicago IL | ³Alverno Laboratories, Hammond IN

BACKGROUND

Matrix-assisted laser desorption/ionization time-of-flight (MALDI-TOF) is a targeted method to detect known genetic variants. Using MALDI-TOF, we initially developed an assay for 119 variants in the *CFTR* gene (Table 1). A preliminary assay was designed, tested for feasibility internally and sent to two volunteer laboratories for evaluation.

MATERIALS & METHODS

Assays for the targeted variants were designed using the online Agena Assay Design Suite (ADS, v2.0, San Diego CA). Primer sequences were aligned in silico against the GRCh38 reference genome. Primer-primer interactions within multiplex reactions were evaluated, and any primers exhibiting predicted interactions were redesigned and re-analyzed. This iterative process was repeated until no predicted interactions remained among primers within the same reaction well. A total of six wells were required to multiplex all variants of interest.

Genomic DNA is subjected to multiplex PCR amplification followed by shrimp alkaline phosphatase (SAP) to remove unincorporated dNTPs. In a subsequent reaction, probes that are adjacent to the variant of interest anneal to the PCR products and undergo single base extension. The detection platform uses matrix-assisted laser desorption/ionization time-of-flight mass spectrometry technology to determine the extension oligonucleotides mass and differentiate variants.

Reference materials were sourced from Coriell Cell Institute for Medical Research (Camden NJ) and Maine Molecular Quality Control, Inc. (Saco ME). DNA extracted from dried blood spots, buccal swabs, saliva, and whole blood was used for the studies. Feasibility studies were performed by multiple operators on samples using 2 different Agena MassARRAY Systems (San Diego CA) on different days, using both the 96- and 384-well format. The following parameters were guard-banded to determine optimum parameters and robustness: annealing temperature (58-62°C), PCR ramp rate (3-6°C/s), PCR reaction mix (0.5-1.5X), unincorporated extension primers (0.5-1.5X), MgCl₂ concentration (4.0±0.5mM) and recommended sample inputs (2-50ng).

Testing was performed by 2 external laboratories to determine accuracy and analytical sensitivity and specificity from orthogonally tested samples (e.g., next generation sequencing and MassARRAY). The external sites used DNA extracted from whole blood using a Roche MagNA Pure (Pleasanton CA) and dried blood spots using a laboratory-developed assay extraction method. Thermal cyclers used included: Applied BioSystems Incorporated (ABI) ProFlex, ABI Verti Pro, and ABI Verti (ThermoFisher, Waltham MA).

Table 1. *CFTR* 100+ Panel Content


Legacy Name	rsID#	HGVS cDNA (NM_00492.4)	HGVS (GRCh38) (NC_000007.14)	ACMG Status
G2X	rs397508740	c.40C>T	g.17480094C>T	Recommended
I506V	rs72734892	c.178G>T	g.17509670G>T	Recommended
P67L	rs368505753	c.200C>T	g.17509696C>T	Recommended
R75X	rs121908749	c.222C>T	g.17509902C>T	Recommended
Q88E	rs79961395	c.254G>A	g.17509123G>A	Recommended
394delTT	rs121908750	c.262_263del	g.17509131_17509132del	Recommended
Q91R	rs121908750	c.271G>A	g.17509140G>A	Recommended
406-1G>A	rs121908792	c.274-1G>A	g.17530889G>A	Recommended
Q96X	rs397508461	c.292C>T	g.17530917C>T	Recommended
Q98H	rs397508464	c.293A>G	g.17530918A>G	Recommended
444delAA	rs121908801	c.313del	g.17530936del	Recommended
D110H	rs113993958	c.328G>C	g.17530953G>C	Recommended
R117H	rs7865421	c.348C>T	g.17530974C>T	Recommended
R117H	rs7865421	c.350G>A	g.17530975G>A	Recommended
G2V1G>T	rs78755691	c.489T>G>T	g.17531150G>T	Recommended
F191V	rs141482808	c.571T>G	g.17534357T>G	Recommended
711+1G>T	rs77188391	c.579+1G>T	g.17534366G>T	Recommended
711+3A>G	rs397508761	c.579+3A>G	g.17534368A>G	Recommended
L208W	rs121908752	c.617T>G	g.17535385T>G	Recommended
L218X	rs397508777	c.653T>A	g.17535321T>A	Recommended
V232D	rs397508783	c.695T>A	g.17535363T>A	Recommended
935delA	rs121908772	c.803del	g.17536607del	Recommended
Q290X	rs397508808	c.886C>T	g.17536672C>T	Recommended
G338K	rs7961340	c.898G>T	g.17536742G>T	Recommended
R338W	rs121909011	c.1000C>T	g.17540230C>T	Recommended
T388I	rs7409459	c.1013C>T	g.17540243C>T	Recommended
I154insTC	rs397508600	c.1021_1022dup	g.17540251_17540252dup	Recommended
I161delC	rs121908774	c.1036del	g.17540273del	Recommended
R347P	rs77932196	c.1040G>C	g.17540270G>C	Recommended
R352Q	rs121908753	c.1055G>A	g.17540285G>A	Recommended
L288insTA	rs121908785	c.1152_1156dup	g.17540264_17540265dup	Recommended
1461insA	rs397508189	c.1322_1330dup	g.17540758_17540761dup	Recommended
A455E	rs74551128	c.1364C>A	g.17548795C>A	Recommended
V456A	rs19322500	c.1367T>C	g.17548798T>C	Recommended
I504delG	rs397508196	c.1373del	g.17548836del	Recommended
I526-1G>A	rs397508209	c.1388-1G>A	g.17549463G>A	Recommended
S466X	rs121908805	c.1397C>A	g.17559468C>A	Recommended
L467P	rs13957311	c.1460T>C	g.17559471T>C	Recommended
I507del	rs121908745	c.1519_1521del	g.17559590_17559592del	Recommended
I508del	rs397508660	c.1522_1523del	g.17559593_17559594del	Recommended
C524X	rs121908754	c.1572C>A	g.17559643C>A	Recommended
L716+1G>A	rs397508230	c.1584+1G>A	g.17559666G>A	Recommended
L717-1G>A	rs76713772	c.1585-1G>A	g.17559778G>A	Recommended
G542X	rs113993959	c.1616G>T	g.17559778G>T	Recommended
S549N	rs121908755	c.1646G>C	g.17559780G>C	Recommended
S549R	rs121909005	c.1648A>C	g.17559782A>C	Recommended
G551S	rs121909013	c.1651G>A	g.17559785G>A	Recommended
G551D	rs75527207	c.1652G>A	g.17559786G>A	Recommended
R552X	rs74507245	c.1653C>T	g.17559787C>T	Recommended
L558S	rs19392204	c.1673T>C	g.17559827T>C	Recommended
A559T	rs75549581	c.1675G>A	g.17559829G>A	Recommended
R567T	rs8005510	c.1679G>C	g.17559833G>C	Recommended
I568del	rs397508623	c.1679-1G>A	g.17559834del	Recommended
I181+1G4A>G	rs397508266	c.1680-888A>G	g.17559847A>G	Recommended
R565G	rs397508267	c.1680A>C	g.17559848A>C	Recommended
A561E	rs121909047	c.1682C>A	g.17559935C>A	Recommended
I182delA	rs19322505	c.1682del	g.17559935del	Recommended
V560D	rs397508276	c.1705T>G	g.17559937T>G	Recommended
I585X	rs397508296	c.1753G>T	g.17559942G>T	Recommended
I898+1G>A	rs121908748	c.1766+1G>A	g.17599480G>A	Recommended
I898+5G>T	rs121908796	c.1766+5G>T	g.17599480G>T	Recommended
R1009T	rs20197862	c.1818G>T	g.17599508G>T	Recommended
G628R	rs397508116	c.1882G>A	g.17599504G>A	Recommended
L284insA	rs121908746	c.2052dup	g.17599219dup	Recommended
L284delA	rs121908746	c.2052del	g.17599219del	Recommended
L290insA	rs121908760	c.2123dup	g.17599229dup	Recommended
L290insA	rs746418935	c.2175dup	g.17599342dup	Recommended
R764X	rs121908810	c.2200C>T	g.17599427C>T	Recommended
R765X	rs374948172	c.2233C>T	g.17599520C>T	Recommended
R792X	rs164584946	c.2233G>T	g.17599520G>T	Recommended
2622+1G>A	rs141158996	c.2499+1G>A	g.17592558G>A	Recommended
L2789+5G>A	rs80224560	c.2657+5G>A	g.17602868G>A	Recommended
Q890X	rs79633941	c.2668C>T	g.17603542C>T	Recommended
I913L	rs149796277	c.2729T>A	g.17602867T>A	Recommended
R552X	rs397508442	c.2738C>T	g.17602868C>T	Recommended
G979D	rs386134230	c.2909G>A	g.17606674G>A	Recommended
L320del-A	rs121908797	c.2988del	g.17606733del	Recommended
L320+1G>A	rs7590551	c.2988+1G>A	g.17606734G>A	Recommended
L319delA	rs121908787	c.302_307del	g.17610997_17610997del	Recommended
L103delN	rs397508498	c.3107C>A	g.17610637C>A	Recommended
L322-2del-A	rs76151804	c.3140-2del	g.17611555del	Recommended
R1066G	rs78194216	c.3196G>T	g.17611637G>T	Recommended
R1064A	rs121908719	c.3197G>A	g.17611638G>A	Recommended
W1089X	rs78802634	c.3266G>A	g.17611707G>A	Recommended
W1098C	rs397508533	c.3294G>C	g.17611735G>C	Recommended
S1118F	rs146521846	c.3333C>T	g.17611794C>T	Recommended
R1158I	rs7969223	c.3472C>T	g.17627325C>T	Recommended
R1162K	rs74647930	c.3484C>T	g.17627325C>T	Recommended
R559delC	rs7898473	c.3528del	g.17627581del	Recommended
W1204X	rs121908765	c.3611G>A	g.17627664G>A	Recommended
L791delC	rs121908811	c.3659del	g.17627712del	Recommended
3849+5G>A	rs19322520	c.3717+5G>A	g.17627756G>A	Recommended
3849+10del-C>T	rs75939782	c.3718-2477C>T	g.17639961C>T	Recommended
3876delA	rs121908784	c.3744del	g.17642464del	Recommended
L1255X	rs76649725	c.3764C>A	g.17642484C>A	Recommended
394delG	rs76384946	c.3888del	g.17642536del	Recommended
W1282X	rs7010888	c.3845G>A	g.17642565G>A	Recommended
4016insT	rs121908808	c.3889dup	g.17642875dup	Recommended
N1303K	rs80034486	c.3909C>G	g.17642877C>G	Recommended
CFTRdel62.3	structural variant	c.653+154-11,3273+1_274-15del	NA	Additional
E20V	rs121908751	c.270A>T	g.17530899G>T	Additional
G178R	rs80282562	c.532G>A	g.17534318G>A	Additional
1078delT	rs75528968	c.948del	g.17540178del	Additional
I38K	rs397508139	c.1007T>A	g.17540237T>A	Additional
Q403X	rs77101217	c.1470C>T	g.17559586C>T	Additional
poly T (S1)	rs1805177	c.1210-7_1210-6del	g.17548634_17548635del	Additional
poly T (O1)	rs1805177	c.1210-7_1210-6dup	g.17548634_17548635dup	Additional
I06V	rs1800991	c.1516A>G	g.17559597A>G	Additional
I07V	rs1801178	c.1519A>G	g.17559598A>G	Additional
F508C	rs74571530	c.1523T>G	g.17559594T>G	Additional
1677delTA	rs121908776	c.1545_1546del	g.17559616_17559617del	Additional
V520F	rs7944904	c.1558G>T	g.17559629G>T	Additional
I181+1C>G	rs397508483	c.1679+1G>C	g.17587834G>C	Additional
G628R-C>C	rs397508116	c.1882G>C	g.17592698G>C	Additional
2055delG>A	rs121908779	c.1923_1931delinsA	g.17592090_17592098delinsA	Additional
L243delT	rs121908812	c.2012del	g.17592179del	Additional
L218del-A	rs121908799	c.2051_2052delinsG	g.17592218_17592219delinsG	Additional
F725K	rs199791061	c.2173G>A	g.17592340G>A	Additional
Y1092X-C>A	rs121908761	c.3276C>A	g.17611717C>A	Additional
Y1092X-C>G	rs121908761	c.3276C>G	g.17611717C>G	Additional
M151K	rs36210737	c.3382T>A	g.17611737T>A	Additional
D1152H	rs75541569	c.3454G>C	g.17614699G>C	Additional
3905insT	rs121908789	c.3773dup	g.17642493dup	Additional

RESULTS

The assay performance during guard-banding was optimal at recommended performance parameters, with suboptimal performance for DNA concentration at 2ng and low MgCl₂ concentration, especially for 394delTT.

The external laboratory results (68 samples) were 100% concordant with previously tested results, for 100% accuracy. Three samples failed testing. Fifty-two chromosomes where the reference was positive tested positive and 78 chromosomes where the reference was negative tested negative. Figure 1 includes the variants included in the orthogonal testing. The analytical sensitivity is 100% (95% CI; 93.1-100%) and the analytical specificity is 100% (95% CI; 95.3-100%).

Figure 1. Variants included in orthogonal testing



ACKNOWLEDGEMENTS

The entire Agena Team: Aaron Florece, Aleksey Nakorchevsky, Amber Lane, Glenn Wise, Maggie Lowe, Heath Metzler, Nichole Rupp, Sharron Ohgi, Shelley Spisak, Wayne Ge, Darryl Irwin, Paul Yates.

DISCLOSURES

VMP, AB, KG, LM, PH, PS, MJ, and GS are employees of Agena Bioscience.

CONCLUSIONS

MALDI-TOF is a targeted method to detect known variants. It has mid-range multiplexing capabilities, allowing the detection of approximately 25 variants/well. As such, we developed a targeted assay to detect over 100 *CFTR* variants. The assay was designed to detect the 100 variants as recommended by ACMG, plus additional variants that were recommended by a voice of the customer survey.

Overall, using MALDI-TOF to detect *CFTR* variants performed well in initial studies. One of the limitations of this study is that not all the variants in the assay could be evaluated in clinical samples or reference materials. Currently there are no commercially available reference materials to verify the accuracy of performance for all the variants in the assay. While a commercial control from Maine Molecular Quality Control, Inc. and Coriell DNA samples from the CDC's-GeT-RM studies were used, synthetic oligos had to be created for many of the variants.

Orthogonal studies performed in the external laboratories were 100% concordant with the previous results. A total of 68 DNA samples were tested. The most common variant in the data set, excluding 5T, was as expected F508del (n=7), followed by R117H (n=4). Note that the samples included in the orthogonal testing were not random but were chosen to include as many variants as possible.

Initial customer feedback indicated that the assay performed well using DNA extracted from a variety of sample types. Customers also suggested adding three variants (I506V, I507V, and F508C) that might interfere with targeted testing, to ensure accurate genotyping. These variants are not pathogenic but will ensure that the assay can accurately detect F508del and I507del when present.

Additional studies are underway both internally and externally to improve and further delineate assay performance.

REFERENCE

JL Deignan, AR Gregg, WW Grody, MH Guo, H Kearney, KG Monaghan, KS Raraigh, J Taylor, CJ Zepeda-Mendoza, C Ziats. Updated recommendations for *CFTR* carrier screening: A position statement of the American College of Medical Genetics and Genomics (ACMG), Genetics in Medicine, 2023 25(8):1098-3600.