

Variant List – UltraSEEK™ EGFR Panel

Gene	CDS	Amino Acid Change	COSMIC ID
EGFR	c.2235_2249del(15)GGAATTAAGAGAAGC	p.E746_A750del	6223
	c.2573T>G	p.L858R	6224
	c.2236_2250del(15)GAATTAAGAGAAGCA	p.E746_A750del	6225
	c.2369C>T	p.T790M	6240
	c.2390G>C	p.C797S	5945664
	c.2389T>A	p.C797S	6493937

For Research Use Only. Not for use in diagnostic procedures.

Agena Bioscience, Inc.
4755 Eastgate Mall
San Diego, CA 92121
Phone: +1.858.882.2800

Order Desk: +1.858.202.9301
Order Desk Fax: +1.858.202.9220
orderdesk@agenabio.com
Web: agenabio.com

US +1.877.4.GENOME
EU +49.40.899676.0
AP +61.7.3088.1600
CN +86.21.6427.0566

