



Variant List – Exome QC Panel

SAMPLE QUALITY			SAMPLE IDENTITY
Target	Amplicon Size	Competitive Template Copy Number	Target
POLR2A	100bp	300	rs1045728
POLR2A	100bp	1000	rs10495563
TBP	100bp	3000	rs1065457
ALB	100bp	10000	rs1127379
IPO8	100bp	30000	rs11998387
POLR2A	200bp	300	rs1200349
POLR2A	200bp	1000	rs12594531
TBP	200bp	3000	rs1344
IPO8	200bp	10000	rs17548783
HBG2	300bp	300	rs2246209
GUSB	300bp	1000	rs2273171
AHSP	300bp	3000	rs2301771
RPL27A	300bp	3000	rs3743165
TPT1	300bp	10000	rs3884596
VIM	300bp	30000	rs4478844
POLR2A	400bp	300	rs586421
RPL27A	400bp	1000	rs6420424
HBS1L	400bp	3000	rs6977125
TPT1	400bp	10000	rs7653897
IPO8	400bp	30000	rs773901
PDCD2	500bp	300	rs9131
TBP	500bp	1000	AHSP
PPIA	500bp	3000	AMEL

For Research Use Only. Not for use in diagnostic procedures.

Agena Bioscience, Inc.
4755 Eastgate Mall
San Diego, CA 92121
Phone: +1.858.882.2800

Order Desk: +1.858.202.9301
Order Desk Fax: +1.858.202.9220
orderdesk@agenabio.com
Web: agenabio.com

US +1.877.4.GENOME
EU +49.40.899676.0
AP +61.7.3088.1600
CN +86.21.6427.0566

