



Variant List – Sample ID Panel

Assay Name	Marker Name	Location
SIDv1_QC01	albumin	chr4
SIDv1_QC02	albumin	chr4
SIDv1_QC03	albumin	chr4
SIDv1_QC04	albumin	chr4
SIDv1_QC05	albumin	chr4
SIDv1_XY01	AMEL_XY	chrX/Y
SIDv1_XY02	ARSD_XY	chrX/Y
SIDv1_XY03	TGIF2L_XY	chrX/Y
SIDv1_SNP01	rs1005533	chr20
SIDv1_SNP02	rs1024116	chr18
SIDv1_SNP03	rs1028528	chr22
SIDv1_SNP04	rs10495407	chr1
SIDv1_SNP05	rs10771010	chr12
SIDv1_SNP06	rs11781516	chr8
SIDv1_SNP07	rs13050660	chr21
SIDv1_SNP08	rs1335873	chr13
SIDv1_SNP09	rs1357617	chr3
SIDv1_SNP10	rs1360288	chr9
SIDv1_SNP11	rs136337	chr22
SIDv1_SNP12	rs1382387	chr16
SIDv1_SNP13	rs1413212	chr1
SIDv1_SNP14	rs1454361	chr14
SIDv1_SNP15	rs1463729	chr9
SIDv1_SNP16	rs1468118	chr17
SIDv1_SNP17	rs1493232	chr18
SIDv1_SNP18	rs1982986	chr1
SIDv1_SNP19	rs1994997	chr12
SIDv1_SNP20	rs2010253	chr17
SIDv1_SNP21	rs2040411	chr22
SIDv1_SNP22	rs2046361	chr4
SIDv1_SNP23	rs2056277	chr8
SIDv1_SNP24	rs2076848	chr11
SIDv1_SNP25	rs214054	chr6
SIDv1_SNP26	rs2247221	chr18
SIDv1_SNP27	rs2518968	chr15
SIDv1_SNP28	rs251934	chr5
SIDv1_SNP29	rs2714854	chr7
SIDv1_SNP30	rs2831700	chr21
SIDv1_SNP31	rs354439	chr13

Variant List – Sample ID Panel (PAGE 2 OF 2)

Assay Name	Marker Name	Location
SIDv1_SNP32	rs3819854	chr3
SIDv1_SNP33	rs717302	chr5
SIDv1_SNP34	rs727811	chr6
SIDv1_SNP35	rs729172	chr16
SIDv1_SNP36	rs740910	chr17
SIDv1_SNP37	rs8037429	chr15
SIDv1_SNP38	rs826472	chr10
SIDv1_SNP39	rs876724	chr2
SIDv1_SNP40	rs891700	chr1
SIDv1_SNP41	rs901398	chr11
SIDv1_SNP42	rs914165	chr21
SIDv1_SNP43	rs9583190	chr13
SIDv1_SNP44	rs964681	chr10

For Research Use Only. Not for use in diagnostic procedures.

2

Agena Bioscience, Inc.
4755 Eastgate Mall
San Diego, CA 92121
Phone: +1.858.882.2800

Order Desk: +1.858.202.9301
Order Desk Fax: +1.858.202.9220
orderdesk@agenabio.com
Web: agenabio.com

US +1.877.4.GENOME
EU +49.40.899676.0
AP +61.7.3088.1600
CN +86.21.6427.0566

