



Performance Verification Kit

For the MassARRAY® System



Verify the MassARRAY® Instrument Performance

The performance verification kit for the MassARRAY system is an optional kit which tests and documents system performance in less than 1.5 hours. It provides an easy method for laboratories to:

- Confirm the MassARRAY and CPM are operating within specifications
- Satisfy quality system and inspection requirements by providing proof of optimal instrument performance
- Facilitate remote customer support


For Research Use Only.
Not for use in diagnostic procedures.

KIT COMPONENTS

	Component	Shipping & storage conditions
SpectroCHiPs	2x384D or 2x96F	Ambient temp. (15°C to 25°C)
3-point calibrant	1x 600uL	Dry ice (-25°C to -10°C)
Performance verification standard	6x 1600uL or 1x 1600uL	Dry ice (-25°C to -10°C)

WORKFLOW

The kit has a streamlined workflow requiring less than 90 minutes to run with 5 minutes of hands-on time.



Total Time to Run: Approx. 1hr 20mins

Hands-on Time: 5 Min

ORDERING INFORMATION

The performance verification kit is available in both 96 and 384 reaction formats. Two kits are included within the Agena Care Essential service plan and twelve kits are included in the Agena Care Plus service plan. Contact customer service to learn more.

Product	Number of Runs	Format	Catalog No.
Performance Verification Kit - CPM (2x96)	2	96 CPM	13338F
Performance Verification Kit - CPM (2x384)	2	384 CPM	13343D

For Research Use Only. Not for use in diagnostic procedures.

Agena Bioscience, Inc.
4755 Eastgate Mall
San Diego, CA 92121
Phone: +1.877.443.6663

Order Desk: +1.858.202.9301
Order Desk Fax: +1.858.202.9220
orderdesk@agenabio.com
Website: www.agenabio.com

Online Support:
<https://support.agenabio.com>

