



Variant List – UltraSEEK® EGFR Panel v2

Gene	Nucleotide Change	Amino Acid Change	COSMIC ID	Genomic Mutation ID
EGFR	c.2235_2249del(15)GGAATTAAGAGAAGC	p.E746_A750del	6223	COSV51765119
	c.2573T>G	p.L858R	6224	COSV51765161
	c.2236_2250del(15)GAATTAAGAGAAGCA	p.E746_A750del	6225	COSV51765066
	c.2369C>T	p.T790M	6240	COSV51765492
	c.2390G>C	p.C797S	5945664	COSV51766509
	c.2389T>A	p.C797S	6493937	COSV51766493

For Research Use Only. Not for use in diagnostic procedures.

Agena Bioscience, Inc.
4755 Eastgate Mall
San Diego, CA 92121
Phone: +1.858.882.2800

Order Desk: +1.858.202.9301
Order Desk Fax: +1.858.202.9220
orderdesk@agenabio.com
Web: www.agenabio.com

Online Support:
<https://support.agenabio.com>

