

Variant List – ClearSEEK™ PIK3CA v1 Panel

Gene	Nucleotide Change	Amino Acid Change	COSMIC ID	Genomic Mutation ID
PIK3CA	c.1258T>C	p.C420R	757	COSV55874020
PIK3CA	c.1624G>A	p.E542K	760	COSV55873227
PIK3CA	c.1624G>C	p.E542Q	17442	COSV55894248
PIK3CA	c.1633G>A	p.E545K	763	COSV55873239
PIK3CA	c.1633G>C	p.E545Q	27133	COSV55878227
PIK3CA	c.1634A>C	p.E545A	12458	COSV55873209
PIK3CA	c.1634A>G	p.E545G	764	COSV55873220
PIK3CA	c.1634A>T	p.E545V	27155	COSV55892885
PIK3CA	c.1635G>C	p.E545D	27374	COSV55875881
PIK3CA	c.1635G>T	p.E545D	765	COSV55874040
PIK3CA	c.1636C>A	p.Q546K	766	COSV55873527
PIK3CA	c.1636C>G	p.Q546E	6147	COSV55882350
PIK3CA	c.1637A>C	p.Q546P	767	COSV55875400
PIK3CA	c.1637A>G	p.Q546R	12459	COSV55876869
PIK3CA	c.1637A>T	p.Q546L	25041	COSV55877455
PIK3CA	c.3139C>A	p.H1047N	5029128	COSV55912376
PIK3CA	c.3139C>T	p.H1047Y	774	COSV55876499
PIK3CA	c.3140A>C	p.H1047P	249874	COSV55888015
PIK3CA	c.3140A>G	p.H1047R	775	COSV55873195
PIK3CA	c.3140A>T	p.H1047L	776	COSV55873401

For Research Use Only. Not for use in diagnostic procedures.

Agena Bioscience, Inc.
4755 Eastgate Mall
San Diego, CA 92121
Phone: +1.858.882.2800

Order Desk: +1.858.202.9301
Order Desk Fax: +1.858.202.9220
orderdesk@agenabio.com
Web: www.agenabio.com

Online Support:
<https://support.agenabio.com>

